

Installation Instructions

There are two main branches of development, follow the link below for the correct installation or upgrade instructions.

- [DcmSuite 1.19.x or DcmMG/DcmRadiologist 1.1.x](#)
- [DcmSuite 1.20.x or DcmMG/DcmRadiologist 1.2.x](#)

Studio Synthesis
Sistemi Avanzati S.r.l.



Via dell'Oreficeria, 32 – 36100 Vicenza – ITALY
tel. (+39) 0444 565499 – fax (+39) 0444 961515

www.studiosynthesis.it

DCMSuite 1.20.x or DcmMG (DcmRadiologist) 1.2.x

Prerequisites

The following installation dependencies are needed in order to perform a successful installation of the DcmSuite/DcmMG programs:

Operating System

Windows XP Home/PRO/Server

Windows XP Home/PRO/Server with Service Pack 2

Windows 2000 PRO/Server with Service Pack 6

Warning: on Windows VISTA we have installation problems which will be fixed as soon as possible.

NB: if you want to make an installation on older systems (at least NT 4.0) you can try by installing IE6 and MDac 2.7 packages from Microsoft.

Firebird SQL Server 2.0.x and Firebird ODBC driver 1.2.x.x

This is the RDBMS we use, an opensource project available for download at <http://www.firebirdsql.org/>. Just follow the standard default installation for both packages without checking any option.

New Installation

Just launch the setup after installing the over mentioned dependencies.

DcmSuite Upgrade from version 1.17.x or above and DcmMG upgrade

The program switched to a new version of the SQL server, faster and more stable. Before upgrading perform the following steps:

1. Uninstall the Firebird ODBC
2. Uninstall the Firebird SQL Server 1.5: this is because two different firebird services cannot live together, anyway the installation program will warn you if another instance of the server is running.
3. Install the Firebird 2.0.x server
4. Install the Firebird 1.2.x ODBC driver
5. Launch the DcmSuite/DcmMG setup program

DcmSuite Upgrade from version 1.16.x or above

Just launch the setup program after installing the over mentioned dependencies.

DCMSuite 1.19.x or DcmMG (DcmRadiologist) 1.1.x

Prerequisites

The following installation dependencies are needed in order to perform a successful installation of the DcmSuite/DcmMG programs:

Operating System

Windows XP Home/PRO/Server

Windows XP Home/PRO/Server with Service Pack 2

Windows 2000 PRO/Server with Service Pack 6

Warning: Windows VISTA doesn't work

NB: if you want to make an installation on older systems (at least NT 4.0) you can try by installing IE6 and MDac 2.7 packages from Microsoft.

New Installation

Just launch the setup.

Upgrade

Just launch the setup.